

Diamond Saw

Diamond Saw

A black Allied Trim Saw 2 is shown from a three-quarter front view. It has a digital display on the left side showing '2500' in green. A transparent safety guard is partially open, revealing the cutting mechanism. The machine is compact and sits on a flat base.

Tool Type: "Cutting Tool"

Location: "Microfluidics Lab"

Supervisor	Tool Lead
David Bothman	
(805) 893-4125	
bothman@cnsi.ucsb.edu	

Description: "Trim Saw 2"

Manufacturer: "Allied"

About

The Trimsaw 2 is a saw designed for rapid sectioning of printed circuit boards, ceramic substrates and other abrasive non metals. Contains a built in reservoir which can be used to feed lubrication directly to the blade.

Safety Concerns

- Keep hands clear of blade when machine is on
- Use only approved blades and do not cut metals with sectioning blades
- Use coolant or lubrication when making cuts
- Use fixtures and guides whenever possible to reduce risk of slipping into the blade

Training Documentation

[trimsaw2SOP](#)

Detailed Specifications

Motor power: .25 HP Speed range: 500-3000 RPM by 100 RPM increments Cutting platform: 6.75" x 12.5" Blade size: 4"-6"

Reference Documentation

- [Allied blade catalog](#)
 - [UKAM Industrial Superhard Tools](#)
 - [Kool Mist MSDS](#)
-

From:

<https://microfluidics.cnsi.ucsb.edu/wiki/> - Innovation Workshop Wiki

Permanent link:

https://microfluidics.cnsi.ucsb.edu/wiki/doku.php?id=diamond_saw&rev=1672885801

Last update: **2023/01/05 02:30**

