# **Thermoline / Lindberg Blue M Oven**



Tool Type: Oven

Location: Innovation Workshop

**Description:** Oven

Manufacturer: Thermo Scientific

Last updated 1/15/23 Haley

### **About**

# **Safety Concerns**

- High Voltage
- High Temperature
- No combustibles or hazardous materials
- Thermally insulating gloves must be worn when working with hot items. They can usually be found on the bench next to the ovens.
- No food may be placed in laboratory ovens.
- Never heat a sealed container.
- Do not operate the oven if it is damaged, or is not working properly. Please place a note on the oven to alert other users, and notify laboratory staff of the problem.
- Double containment is required for heating containers of liquid to prevent leaked material from contaminating the oven. This is especially true when heating PDMS molds.

## **Operating Procedures**

1. Enter the set point temperature (see instructions for the different controllers below)

#### Last update: 2023/03/15 22:30

- 2. Place your items in the oven. Remember to place liquid materials on a baking sheet or other secondary container.
- 3. Set the countdown timer for the length of your run. If the timer reads zero time remaining another user who needs to use the oven may remove your items.
- 4. Remove your items from the oven and close the door.
- 5. The oven is left powered on.

### Changing the oven set point

- Use the ↑↓ to adjust the temperature
- 2. Press the SET/ENT button

### Setting the countdown timers

- 1. Make sure that the slide switch on the right is on countdown
- 2. Enter the run time (HMS Buttons)
- 3. Press Start/Stop

# **Detailed Specifications**

### **Reference Documentation**

Laboratory Oven SOP

From:

https://microfluidics.cnsi.ucsb.edu/wiki/ - Innovation Workshop Wiki

Permanent link:

https://microfluidics.cnsi.ucsb.edu/wiki/doku.php?id=lindberg-oven

Last update: 2023/03/15 22:30

