

Thermoline / Lindberg Blue M Oven

Thermoline / Lindberg Oven



Tool Type: Oven

Location: Innovation Workshop

Description: Oven

Manufacturer: Lindberg/Blue M

About

Safety Concerns

- High Voltage
- High Temperature
- No combustibles or hazardous materials
- Thermally insulating gloves must be worn when working with hot items. They can usually be found on the bench next to the ovens.
- No food may be placed in laboratory ovens.
- Never heat a sealed container.
- Do not operate the oven if it is damaged, or is not working properly. Please place a note on the oven to alert other users, and notify laboratory staff of the problem.
- **Double containment is required for heating containers of liquid to prevent leaked material from contaminating the oven. This is especially true when heating PDMS molds.**

Operating Procedures

Detailed Specifications

Reference Documentation

[Laboratory Oven SOP](#)

[lindbergovenmanual.pdf](#)

From:

<https://microfluidics.cnsi.ucsb.edu/wiki/> - Innovation Workshop Wiki

Permanent link:

<https://microfluidics.cnsi.ucsb.edu/wiki/doku.php?id=lindberg-oven&rev=1672955542>

Last update: **2023/01/05 21:52**

