

Programming Handout - Thinky ARE-310 Mixer

Last Edited: Bassari 3/3

The thinky mixer is one of the PDMS tools in the Microfluidics lab. This machine is used for mixing and de-aerating the elastomers that are used in fabrication of Microfluidic devices. The mixing time and revolution speed of the device can be programmed using the machine's controller interface.



I) Operation Mode: The device can operate in standard (STD) and stepwise (STEP) modes. Stepwise mode allows access to multiple steps of operation and enables the user to vary the revolution speed, while standard mode only allows for variation of mixing time within one step. You can switch between these modes using the **MODE** button.

The STD recipe works well for mixing Sylgard 184 resin. It should be programmed for 30 seconds of mixing at 2000 RPM and 30 seconds of defoaming at 2200 rpm. **Please do not change this recipe.** If you would like to use different settings please use the STEP mode.

II) Mixing and Defoaming: Each step of operation can have one or two sub-steps - mixing and, possibly, defoaming. The MIX and DEFOAM LEDs highlight the sub-step for which the time and revolution speed are being displayed.

III) Changing Operating Conditions:

To change the time and revolution speed for your process,

- A) Navigate to the step that you want to alternate using the **STEP** button. The current step number is displayed on the step LED's. Usually, processes only use the first step (mixing and defoaming constitute one step).
- B) Use the **SETUP** button to cycle through the editable fields for that step - these are, in order of appearance, time of mixing, speed of mixing, display of mixing settings, time of defoaming, speed of defoaming, and display of defoaming settings followed by display of mixing settings and a return to the beginning of the cycle.
- C) At each editable field, alter the numerical value using the **UP** and **DOWN** arrows, then press **SETUP** to move to the next field.
- D) After making changes to the mixing settings, hold the **MEMORY** button for 1+ seconds to save your changes to the step.
- E) After making changes to the defoaming settings, hold the **MEMORY** button for 1+ seconds to save your changes to the step.

IV) Altering the Memory Slot:

By default, the above operations change the settings of the program saved in the first memory slot, as displayed by the memory LED's. There are up to 6 programs that can be saved on the machine. To use a different slot, press the **MEMORY** button.