

# Lindberg Blue M Convection Oven

## Convection Oven



**Tool Type:** Oven

**Location:** Microfluidics Lab

**Description:** Mechanical Convection Oven

**Manufacturer:** Lindberg / Blue M

## About

## Safety Concerns

High Voltage, High Temperature, no Combustibles or Hazardous Materials

## Operating Procedures

1. Enter the set point temperature (see instructions for the different controllers below)
2. Place your items in the oven. Remember to place liquid materials on a baking sheet or other secondary container.
3. Set the countdown timer for the length of your run. If the timer reads zero time remaining another user who needs to use the oven may remove your items.
4. Remove your items from the oven and close the door.
5. The oven is left powered on.

## Detailed Specifications

## Reference Documentation

[Laboratory Oven SOP](#)

Operations Manual

---

From:  
<https://microfluidics.cnsi.ucsb.edu/wiki/> - **Innovation Workshop Wiki**

Permanent link:  
<https://microfluidics.cnsi.ucsb.edu/wiki/doku.php?id=lindberg-convection-oven&rev=1672952693>

Last update: **2023/01/05 21:04**

